EmbryoGlue®

US

EN: Indication for use

Medium for embryo transfer.

Product Description

EmbryoGlue® is a bicarbonate buffered medium containing recombinant human albumin, hyaluronan and gentamicin as an antibacterial agent.

For use after equilibration at +37°C and 6 % CO₂ atmosphere.

Storage instructions and stability

Store dark at +2 to +8 °C.

EmbryoGlue® is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

For Non-EU: Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

EU: Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record the opening date on the bottle. Discard excess media no later than two weeks after first opening.

Directions for use

Embryo transfer

Add approximately 1 mL of EmbryoGlue® to the well of a rinsed 1-well dish. Add approximately 2 mL of EmbryoGlue® to the moat of the 1-well dish. Equilibrate the dish in +37°C and 6 % CO₂ for 4–24 hours. Equilibrate the embryos to be transferred in the well containing EmbryoGlue® for a minimum of 10 minutes in a 6% CO₂ environment prior to transfer. Rinse the 1 mL non-toxic syringe by drawing up and then out media from the moat several times until no air bubbles are observed in the syringe. Draw up approximately 0.5 mL of the medium from the moat. Firmly attach the transfer catheter to the pre-rinsed 1 mL non-toxic syringe. Flush approximately 0.5–1.0 mL of equilibrated transfer medium from the moat through and out of the catheter. After rinsing, draw approximately 0.1 mL of additional EmbryoGlue® into the moat until approximately 20 µL is left in the syringe. Under microscopic control, gently load the embryos into the catheter in approximately 5–10 µL of additional EmbryoGlue® followed by a small amount of air. (The small pocket of air at the tip allows better imaging of the tip for ultrasound guided transfers). For the embryo transfer, pass the tip of the catheter into the uterus approximately 1 cm from the top of the cavity and expel the embryos in a total volume of approximately 25–30 µL of medium. Slowly withdraw the catheter while maintaining steady pressure on the plunger of the syringe. Make a final microscopic examination of the catheter.

Sterile filtered	SAL 10 ³
Mouse embryo assay (1-cell) [% expanded blastocyst within 96 hours]	≥ 80
Bacterial endotoxins (LAL assay) [EU/mL]	< 0.25

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use EmbryoGlue® if it appears cloudy.

EmbryoGlue® contains recombinant human albumin (produced in yeast), hyaluronan and gentamicin. Do not use in patients with known hypersensitivity to any of the components. To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife’s IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

Description of ISO Symbols

STERILE	Sterilized using aseptic processing techniques
Temperature limit	
Use by - see label.	
Caution: Consult accompanying documents.	
Catalog number	
Batch code	
CE Mark (Conformité Européen)	

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Directions for use

Embryo transfer

Add approximately 1 mL of EmbryoGlue® to the well of a rinsed 1-well dish. Add approximately 2 mL of EmbryoGlue® to the moat of the 1-well dish. Equilibrate the dish in +37°C and 6 % CO₂ for 4–24 hours. Equilibrate the embryos to be transferred in the well containing EmbryoGlue® for a minimum of 10 minutes in a 6% CO₂ environment prior to transfer. Rinse the 1 mL non-toxic syringe by drawing up and then out media from the moat several times until no air bubbles are observed

in the syringe. Draw up approximately 0.5 mL of the medium from the moat. Firmly attach the transfer catheter to the pre-rinsed 1 mL non-toxic syringe. Flush approximately 0.5–1.0 mL of equilibrated transfer medium from the moat through and out of the catheter. After rinsing, draw approximately 0.1 mL of EmbryoGlue® from the center well and expel into the moat until approximately 20 µL is left in the syringe.

Under microscopic control, gently load the embryos into the catheter in approximately 5–10 µL of additional EmbryoGlue® followed by a small amount of air. (The small pocket of air at the tip allows better imaging of the tip for ultrasound guided transfers). For the embryo transfer, pass the tip of the catheter into the uterus approximately 1 cm from the top of the cavity and expel the embryos in a total volume of approximately 25–30 µL of medium. Slowly withdraw the catheter while maintaining steady pressure on the plunger of the syringe. Make a final microscopic examination of the catheter.

Mouse embryo assay (1-cell) [% expanded blastocyst within 96 hours]	≥ 80
Bacterial endotoxins (LAL assay) [EU/mL]	< 0.25

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